ZSTEK ZR5HB v1.10 Ultra-HD Video Imaging System Data Manual

1 Characteristic

- 3-channel 4K60 (3840x2160) coaxial camera input
- 4 HDMI outputs, supports up to 4K60 (3840x2160) resolution
- 5 USB3.0 ports: Support USB storage devices, USB mice, and other peripherals
- 1 Gigabit Ethernet
- ADC input: supports up to 16 key extensions
- 4 RS232 ports: Enables light source (MLS0x) synchronization and serial port screen menu operations
- A clock battery interface to keep the real-time clock running after a power failure
- Power input: 12V 3A(8.5V~13V)
- Size: 154.9mmx110mm
- Multiple metering modes: average, center, and peak
- One-click white balance adjustment with built-in color styles and adjustable color parameters
- Prevents color spillage, aiding in identifying surgical details obscured by blood-red during bleeding procedures.
- Support horizontal and vertical mirroring
- Image freeze, video recording, and snapshots
- Dimming level adjustable, smoke elimination
- Store multiple scenario parameters for easy direct access in different scenarios
- Customize shortcut keys
- The system menu can be hidden and authorized for login, providing a set of serial port commands for system control, enabling users to develop customized serial port screen menu applications.
- Supports web applications and provides SDK packages
- Supports DICOM protocol

2 Summary

The ZR5HB_v1.10 motherboard is a dual-system platform that supports two camera configurations: either a 4K 3D white light system or a 4K 2D fluorescent system, or three cameras with software-switched modes for both 4K 3D white light and 4K 2D fluorescent outputs. Featuring the latest fifth-generation image processing system with advanced ISP capabilities, it delivers real-time low-latency image processing for up to two 4K60 streams. The system incorporates dual-core NNIE processors with 4TOPS computing power, providing essential hardware support for future intelligent image analysis applications. This design ensures sustainable upgrades and facilitates distributed application deployment planning.

For close-up applications such as endoscopy, surgical microscopy, and industrial inspection, the system has been optimized for white balance locking, automatic metering, bright light suppression, smoke elimination, and red coverage in large bleeding environments, ensuring maximum recognition of tissue or object details in the field of view.

The TX22 and TX21A series adapter boards enable integration with OV series electronic microscopes, creating portable hard-tube and electronic microscope applications.


3 Main Application Scenarios

Medical endoscopes Surgical robot navigation Surgical microscopes





4 Apply Block Diagram





Interface Definition

Note: The " \triangle " marks the first pin (Pin1).

name	item	Pin definition	direction for use		
USB 3.0	J22	Type A Dual USB 3.0 Port USB 3.0 port	Supports USB storage device recording and snapshot capture, and USB mouse and keyboard operation		
USB 3.0	J21	Type A Dual USB 3.0 Port USB 3.0 port	Supports USB storage device recording and snapshot capture, and USB mouse and keyboard operation		
Net	J20	1000MHz Ethernet: RJ45	Provide wired Ethernet connection applications		
HDMI	J19	HDMI Type A HDMI OUT	HDMI 1.4: Supports up to 1080P60 resolution output		
HDMI	J18	HDMI Type A HDMI OUT	HDMI 2.0: Supports up to 4K 60Hz resolution output		
HDMI	J17	HDMI Type A HDMI OUT	HDMI 2.0: Supports up to 4K 60Hz resolution output		
Audio	Ј8	Wafer: TH1.25 * 7 PIN1: AC_INR/N PIN2: AC_INL/P PIN3: GND PIN4: Audio_OUT_L PIN5: Audio_OUT_R PIN6: GND PIN7:	Microphone input and audio output		
CMOS4	J24	Wafer: TH1.25 * 3 PIN1: RS232_TXD4 PIN2: RS232_RXD4 PIN3: GND	The system provides a single RS232 standard level interface.		
HDMI	Ј3	HDMI Type A HDMI OUT	HDMI 2.0: Supports up to 4K 60Hz resolution output		
Clock_VCC	J1	Wafer: TH1.25 * 2 PIN1: 3.3V PIN2: GND	The external power supply provides real-time clock power with a maximum voltage of 3.3V.		
RF junctor	J4	Chimeric height 2.5mm 2.5H board end	Bluetooth antenna base		

ZSTEK ZR5HB_v1.10 Ultra-HD Video Imaging System Data Manual

CMOS	J23	TH1.25 * 12H PIN1: RS232_TXD0 PIN2: RS232_RXD0 PIN3: GND PIN4: RS232_TXD1 PIN5: RS232_RXD1 PIN6: GND PIN7: RS232_TXD2 PIN8: RS232_TXD2 PIN8: RS232_RXD2 PIN9: GND PIN10: LIGHT_ADJ PIN11: GND PIN12: 5V/12V	The system provides three RS232 standard level interfaces. • Designed to adjust the parameters of the MLS0xx series light source, it connects to the constant current board (EP1) via a serial port. • Connect the serial port touch screen and access the touch screen operating system menu.
CSI3	J5 J25 J6	J5: MMCX CS5 IN J25: MMCX CS6 IN J6: Wafer: TH1.25 * 3 PIN1: CAM 5V PIN2: GND PIN3: key signal in	Supports dedicated cameras such as QS81/QS82 series
CSI2	J7 J9 J10	J7: MMCX CS3 IN J9: MMCX CS4 IN J10: Wafer: TH1.25 * 3 PIN1: CAM 5V PIN2: GND PIN3: key signal in	Supports dedicated cameras such as QS81/QS82 series
CSI1	J11 J12 J13	J11: MMCX CS2 IN J12: MMCX CS1 IN J13: Wafer: TH1.25 * 3 PIN1: CAM 5V PIN2: GND	Supports dedicated cameras such as QS81/QS82 series

ZSTEK ZR5HB_v1.10 Ultra-HD Video Imaging System Data Manual

		PIN3: key signal in	
		FPC socket: FH41-20S-0.5SH	
		PIN1: ADC	
		PIN2: GND	
		PIN3: USB31_DM	
		PIN4: USB31_DP	
		PIN5: GND	
		PIN6: 5V0_USB3_1	
		PIN7: 5V0_USB3_1	
		PIN8: USB31_RXP	
		PIN9: USB31_RXM	
USB Ext	J14	PIN10: GND	One USB3.0 expansion port and one ADC input
		PIN11: USB31_TXP	
		PIN12: USB31_TXM	
		PIN13: GND	
		PIN14: SYS_PWRON#	
		PIN15: GPIO_UART_TXD	
		PIN16: GPIO_UART_RXD	
		PIN17: 5V	
		PIN18: GND	
		PIN19: 3.3V	
		PIN20: 3.3V	
		3.96 terminal block	To ensure stable system operation, the 12V supply
DC IN	J15	PIN1: GND	must provide 3A drive capability with ripple not
		PIN2: DC12V IN	exceeding 60mVp-p.

6 Electrical Specifications

1) work environment

4				
parameter	minimum	typical case	maximum	explain
input voltage (V)	7.5	12	13.5	
working current		3000		
(mA)				
ambient	5	25	45	
temperature (°C)				
humidity (%rh)	5	55	85	
atmos (kPa)	50	101	106	




Electromagnetic compatibility (YY9706.102-2021:6.1.1 Electromagnetic Interference: Conductive and 2) Radiated Emissions; 6.2.x Antenna Interference)


Test conditions: AC 220V/50Hz, ZR5H (ZS04E metal housing)

electromagnetic distrubance:

Conduction emission: Meets the GB 4824-2019 power terminal electromagnetic interference (EMI) voltage limit (Group 1, Class B) requirements



Radiation emission: Meets the radiation disturbance (Group 1, Class A) limit requirements specified in GB 4824-2019





ZR5HB_v1.10 Ultra-HD Video Imaging System Data Manual

B. Robustness: Meets the performance criteria specified in YY 9706.102-2021

			1
YY 9706.102-2021	project	grade	bear
			fruit
(22	electrostatic discharge	Contact discharge (±6 kV)	
6.2.2		Air discharge (Air): ±8KV	
6.2.3	Radio frequency		
	electromagnetic field	3V/m, 80%AM (1kHz)	qualified
	radiation immunity		
		At the AC power supply port:	
6.2.4	Fast transient pulse cluster	Test voltage peak ±2 kV, repetition frequency 5 kHz or 100 kHz,	qualified
		5/50 ns Tr/Td waveform	
		At the AC power supply port:	
		Line-to-line: ±1 kV peak voltage, open-circuit voltage waveform	
6.2.5	surge	1.2/50µs	qualified
		Line-to-ground: ±2 kV peak voltage, open-circuit voltage	
		waveform 1.2/50µs	
6.2.6	RF field induced conducted	237000/ 434 (1111)	
	immunity	3V,80%AM (1kHz)	qualified
	The voltage on the power	1. Test voltage <5%Ut, lasting 0.5 cycles	
	supply input line is	2. Test voltage at 40% of Ut for 5 cycles	
6.2.7	temporarily	, i	qualified
	Down, short interruptions,	3. Test voltage at 70% of Ut for 25 cycles	
	and voltage changes	4. Test voltage <5%Ut, for 250 cycles	
6.2.8.1	Power frequency magnetic		
	field	3A/m	qualified

7 Related components: SS82, QS82T, MLS0x, TX23